


# Standard Apron Units and Utility Tables

These apron units provide knee space under laboratory work areas so that the laboratory worker may work from the sitting position. The aprons also provide exposed storage areas for large equipment or materials. The standard units are made to the standard depth, heights and lengths of the base units so they can be easily integrated into the laboratory layout.


All aprons are 22" d x 4 3/8" high  
 Note: Suffix "-28" for 28" deep aprons




Type A unit



Type A-1  
one drawer unit






Type A-2  
two drawer unit



Type A-3  
three drawer unit

UNIT	CATALOG NUMBER	OVERALL APRON LENGTH	SIZE OF DRAWER(S)	UNIT	CATALOG NUMBER	OVERALL APRON LENGTH	SIZE OF DRAWER(S)
A	A-324	24 1/8		A-1	A-324-1	24 1/8	F
	A-330	30			A-330-1	30	F
	A-335	35			A-335-1	35	E
	A-342	42			A-342-1	42	E
	A-347	47			A-347-1	47	E
	A-352	52			A-352-1	52	E
	A-358	58			A-358-1	58	E
	A-370	70			A-370-1	70	E
A-2	A-342-2	42	F	A-3	A-342-3	42	F
	A-347-2	47	F		A-347-3	47	F
	A-352-2	52	F		A-352-3	52	F
	A-358-2	58	E		A-358-3	58	F
	A-370-1	70	E		A-370-3	70	F

All length sizes in inches, Drawer: Sizes (I.D.) .  
 F= 13 5/8"W x 2 1/8"H x 19 1/8"D,  
 E= 13 5/8"W x 2 1/8"H x 19 1/8"D



31" High Number	37" High Number	L	W	31" High Number	37" High Number	L	W
750-301	750-351	30	24	750-316	750-366	60	24
750-302	750-352	30	30	750-317	750-367	60	30
750-303	750-353	36	24	750-318	750-368	60	42
750-304	750-354	36	30	750-319	750-369	60	48
750-305	750-355	42	24	750-320	750-370	66	24
750-306	750-356	42	30	750-321	750-371	66	30
750-307	750-357	42	42	750-322	750-372	66	42
750-308	750-358	48	24	750-323	750-373	66	48
750-309	750-359	48	30	750-324	750-374	72	24
750-310	750-360	48	42	750-325	750-375	72	30
750-311	750-361	48	48	750-326	750-376	72	42
750-312	750-362	54	24	750-327	750-377	72	48
750-313	750-363	54	30				
750-314	750-364	54	42				
750-315	750-365	54	48				

31" High Number	37" High Number	L	W	No. of Drawers
750-328	750-378	30	24	1
750-329	750-379	30	30	1
750-330	750-380	36	24	1
750-331	750-381	36	30	1
750-332	750-382	42	24	1
750-333	750-383	42	30	1
750-334	750-384	48	24	2
750-335	750-385	48	30	2
750-336	750-386	54	24	2
750-337	750-387	54	30	2
750-338	750-388	60	24	2
750-339	750-389	60	30	2
750-340	750-390	66	24	3
750-341	750-391	66	30	3
750-342	750-392	72	24	3
750-343	750-393	72	30	3

The catalog number includes a 1" thick Phenolic Core counter top, other tops available when specified. The tables numbered 750-301 through 750-327 and 750-351 through 750-377 can be made of a all-welded 2" tube construction. To indicated all welded construction add "AW" at end of number, i.e.: 750-301-AW. Vibration pads between work surface and frame can be provided when specified. To indicated balance table add "BT" to end of part number, i.e.: 750-301-BT.