


Electric Drying Cabinets

Electric Drying Cabinets are of double wall construction and fully insulated. Five expanded adjustable metal shelves are furnished in addition to the expanded metal shelf over the electric strip heaters in the bottom of the cabinet. A thermostatic control, 0-250° F., 25 Amp, 220 volts AC for the strip heaters is included. DE-472 has a power consumption of 3000 watts, single phase 220 volts and the DE-372 has a power consumption of 3500 watts, single phase, 220 volts.

CATALOG NO.	LENGTH	DEPTH	HEIGHT
DE-372	35"	22"	84"
DE-472	47"	22"	84"




Acid Storage Cabinets

Acid Storage Cabinets provide a variety of storage space for reagents and solvents. They blend with other laboratory equipment and furniture in the laboratory. These cabinets have 5 adjustable shelves and will accommodate many different sizes of containers.

CATALOG NO.	LENGTH	DEPTH	HEIGHT
R8-121	35"	16"	84"
R8-122	47"	16"	84"

*Special depth of 22" is available, when specified. Indicate with "22" suffix, i.e. "R8-123-22"

Special Purpose Storage



Wardrobe Cabinets

Wardrobe storage cabinets are 22" deep and 84" high. These cabinets are provided with a horizontal bar for hanging clothes and a fixed hat shelf on one half of the cabinet. The other half comes with 5 adjustable shelves. (FD-272W 24-1/8" has only hat shelf and clothes rod.)

CATALOG NO.	LENGTH	DEPTH	HEIGHT
FD-272W	24-1/8"	22"	84"
FD-372W	35"	22"	84"
FD-472W	47"	22"	84"



Standard Storage Assemblies

The R-6 Series of storage cabinet assemblies combines the features of wall cases with sectional unit storage and a convenient work ledge. Upper section contains 3 adjustable shelves. Base cabinet contains 1 adjustable shelf. Locks can be furnished, when specified. Work ledge has a 1 inch overhang on 3 sides.

R-6 Type 1-Framed Sliding Glass Doors Above
Steel Hinged Doors with 2 Drawers Below

CATALOG NO.	LENGTH	DEPTH	HEIGHT
R7-121	37"	23"	84"
R7-122	49"	23"	84"

* Other cabinet combinations available upon request.